

First time using Tasso?



Watch our instructional video:
www.tassoinc.com/tassoplus-video



or scan the QR code
to access.

Contact us for help:
1-800-257-2370

Additional Supplies



Timer



Mirror



Pen



Razor



1 Prepare

Step 1



Wash hands before use. Healthcare workers must wear gloves during collection.

Step 2



Identify collection site – for best results remove any hair.

Step 3



Collect all needed items. Remove device from tray.

Step 4



Twist and remove cap from tube.

Step 5



With fill lines facing out and label facing back, press tube into device until snug.

Step 6



Follow heat pack instructions to activate. If pack does not heat up, discard and use hand instead.

Step 7



Increase blood flow at collection site.
Rub heat pack or hand up and down user's shoulder for 2 minutes to warm and increase blood flow.

Step 8



Clean area with alcohol pad and allow to dry.

Step 9



Remove clear cover from red button.

2 Collect

Step 10



Peel tab away from device and stick device to upper arm. Hang arm straight down at side.

Step 11



Press button all the way down ONCE and release. Set timer for 5 minutes. Watch tube fill using mirror as needed.

Step 12



You may not feel the device. Leave on 5 minutes or less if blood at fill line.

Blood may not appear for first 1-2 minutes. Allow tube to fill for 5 minutes. Remove device if blood reaches top line before time is up. Do not fill past top line!

Step 13



Remove device by peeling from one side. Apply supplied bandage vertically and add slight pressure.

Step 14



Remove tube from device with a slight twist and pull down.

Step 15



Snap cap fully onto tube.

Step 16



Gently invert tube 10 times to mix the blood. Tap or flick tube to get blood moving. Any gel in tube will not mix with blood.

3 Ship

Step 17



Write user's date of birth (MM/DD/YYYY) on tube label in the space above the barcode.

Step 18



Place tube in clear return bag and press to seal. Do not remove pad.

Step 19



Record time and date on return bag label.

Step 20



Follow provided instructions to return user's sample to the lab.

Step 21



Follow local guidelines for discarding all materials.

Step 22



Wash hands after use.

